Dipstick BENCE JONES Free K/λ Test

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Dipstick BENCE JONES Free K/λ Test

Competitive immunochromatographic test for the detection of free light chains K and λ in urine

20 det Ref GD68000 00

IVD

LOT on labels

2 - 8°C

SUMMARY AND EXPLANATION

The presence of light-chains specific free K/λ in urine may be an indication of a mono or polyclonal gammopathy. This pathology, called Bences Jones Proteinuria, is characterized by the indiscriminate proliferation of B-cell lymphocytes which are detectable in the urine as free K & λ proteins. The pathology is associated with Multiple Myeloma, Waldenström's Macroglobulinemia, monoclonal light chain amyloidosis and light chain deposition disease.

PRINCIPLE OF METHOD

The MINIAS GLOBE DIAGNOSTICS Dipstick Bence Jones Free K/λ kit is a competitive, immuno-chromatographic assay, based on lateral flow nitrocellulose membrane technology. The assay has 3 reaction zones without prozone effect:

1. The first zone (closest to dipstick tab) represents the procedural control (CP).
2. The second zone (middle) line striped with free λ antigens.
3. The third zone (lowest) line striped with free K antigens.

During the incubation of the dipstick device in unconcentrated, diluted urine, free light chains bind to the colloidal gold conjugated antibodies present in the nitrocellulose membrane, resulting in the formation of antigen-conjugate complexes in the presence of free light chain proteins. During the lateral flow the sample starts to migrate along the three test zone and the free light chain proteins, if present in the sample, will be bound to the specific conjugated antibody. The sample continues to migrate to the control zone which is specific against the capture conjugate to indicate proper test function.

CONTENTS OF KIT

Dipstick - containing 2 desiccants

Quantity: 1x20 det Ref CD20 Ref US20MI

Sample Diluent – vial containing pH 9.2 buffered solution 10mM TRIS, with surfactants and preservative Sodium Azide 0,09%.

Volume: 3 ml Ref DCBJ03 Ref US20MI

Product Instructions – (Ref INSBJMI) - Present document

MATERIALS NECESSARY BUT NOT PROVIDED

- Clean test tubes and tube rack.
- Precision pipettes with disposable tips for dispensing 100 and 200µl.

CONSERVATION AND STABILITY OF REAGENTS

- The kit must be stored at 2 - 8°C and use not after the expiry date printed on the labels.
- Rapidly remove the required dipsticks from their wrapping and avoid extended exposure to the ambient relative humidity.

CONSERVATION AND STABILITY OF SAMPLES

- Use the first-morning- urines or within 24 hours collection.
- As the presence of endogenous/exogenous protease in the sample may lead to the denaturation and fragmentation of the Ig it is recommended to add an anti-protease inhibitor or Sodium Azide 1%, at the beginning of the urine collection.
- Samples must be stored at 2 - 8°C for 2 days or at - 20°C for 1 month, thaw only once.
- Turbid urine samples need sedimentation.

SAFETY PRECAUTIONS

- The reagents are not considered dangerous according to the 67/548/EEC, 88/379/EEC and 76/769/CEE directives about classification, packaging and labelling of dangerous substances. However, they should be handled with caution, avoiding swallowing and contact with skin, eyes and mucous membranes.
- All reagents contained in the kit are for *in vitro* diagnostic use.
- At all times follow Good Laboratory Practice (GLP) guidelines.
- Handle all test specimens with utmost care to avoid infection with potentially hazardous fluids. Avoid skin contact by wearing protective disposable gloves.
- Sodium Azide (NaN₃) contained in some of the kit reagents may react with lead or copper plumbing and form highly explosive metal azides. When disposing reagents flush with large amounts of water to avoid a build-up.
- Dispose of all used materials into appropriate bio-hazard waste containers.

TECHNICAL PRECAUTIONS

- Follow carefully instructions for use described in the test procedure.
- Do not use reagents beyond the expiry date printed on each component label and do not mix reagents from different kit lots.
- Avoid all cross-contamination of samples by using disposable pipette tips.
- Prior to their use, let all reagents equilibrate 30 minutes at ambient temperature.


INTENDED END USER

This diagnostic device is to be used exclusively by laboratory diagnostic professionals (graduated or qualified technicians in Microbiology and Immunochemistry).

PROCEDURE

1. Select urine samples to be tested and equilibrate at room temperature.
2. Prepare equivalent number of sample tubes and mark them with patient ID's.
3. Pipette into each sample tube 100 µl of Sample Diluent.
4. Pipette into each sample tube 100 µl of urine.
5. Insert one dipstick into each sample tube keeping it with fingers by the MGD logo side (see image below).
6. Incubate 15 minutes (±1 min.) at room temperature.
7. Read and match results comparing them with *Reading-ruler* included in the kit.


Notes – Important! - Do not read test results over 20 minute, beyond this time results are invalid; - Do not touch membrane with finger, handle the dipstick by MINIAS K/λ logo side;

PIPETTE 100 µl URINE SAMPLE AND 100 µl SAMPLE DILUENT

INSERT THE DIPSTICK INTO SAMPLE TUBE

RESULTS INTERPRETATION

Ref	Rev	Dat	Pag
GD680000	01	2011/12/01	3 / 4

ACCEPTABILITY CRITERIA


The purple line closest to the dipstick tab, representing the procedural control, must be visible to consider the assay result valid.

INTERPRETATION OF RESULTS

Use the Reading-ruler included in the kit.

- **Positive K:** Presence of purple line in the second reaction zone of membrane (plus procedural control line) (1).
- **Positive λ :** Presence of purple line at the lowest reaction zone of membrane (plus procedural control line) (2).
- **Positive K e λ :** Presence of purple line at the first reaction zone of membrane only (procedural control line) (3).
- **Negative K e λ :** Presence of purple line in second and third reaction zone (plus procedural control line) (4).
- **Invalid:** Absence of three purple lines.

ATTENTION: a different colour intensity of the procedural control line respect to the reaction zones does not influence the final result and its interpretation.

POS K 1 POS λ 2 POS K λ 3 NEG 4 INVALID 5

TROUBLESHOOTING

1. **Sample migration is not homogeneous:**
 - adsorbent pad and white adhesive at dipstick base-hedge are not fully fixed onto the membrane: repeat the assay with a new dipstick
 - adsorbent pad and coloured/serigraphated adhesive at the dipstick top-hedge are not fully fixed onto the membrane: repeat the assay with a new dipstick
 - the membrane is damaged: repeat the assay with a new dipstick
 - sample volume is not adequate: repeat the assay using the correct amount of sample and diluent
 - turbid urine sample: centrifugate urine and repeat assay
2. **Bands do not appear - procedural control do not appear**
 - Reagents and/or sample temperature is too low: allow all reagents and sample to reach room temperature
 - adsorbent pad and white adhesive at dipstick base-hedge are not fully fixed onto the membrane: repeat the assay with a new dipstick
 - the membrane is damaged: repeat the assay with a new dipstick
 - sample volume is not adequate: repeat the assay using the correct amount of sample and diluent
3. **Bands shape is not well defined:**
 - storage condition of reagents has been not correct: store at 2-8° C;
 - reagents has been used beyond expire date: repeat assay using a new kit
 - see point 1.
4. **Bands appear after 20'**
 - Results obtained beyond 20 minutes are not valid

TEST PERFORMANCE

Sensitivity – The sensitivity of the MINIAS GLOBE DIAGNOSTICS Dipstick Bence Jones assay was assessed testing positive K and λ serially

diluted urine samples, by nephelometry: the sensitivity of the test is 30 mg/l for K and λ free light chains.

Specificity – A study with 10 negative, 10 K and 10 λ positive samples showed 100 % specificity.

Reproducibility – A series of negative and positive K and λ samples were repeatedly tested and acceptable inter-and intra assay CVs were obtained.

LIMITATION OF PROCEDURE

It has been demonstrated that high levels of albumin, glucose and organic acids don't interfere with the assay.

Also pH variations don't interfere with the final test result.

WASTE DISPOSAL

Refer to local law.

BIBLIOGRAFIA/REFERENCES

1. Girkontatie I et al. (1996) A rapid ELISA test for the detection of human Paraproteins, Eur. J. Clin. Chem. Clin. Biochemm 1996 Apr; 34 (4):349-353.
2. Guinan JE et al. (1986) Detection and typing of serum paraproteins with the Quantimetric kappa: lambda ratio test. Clin. Chem. 1986 Oct; 32 (10) 1981-1982.
3. Rowe DS (1968) Quantitative estimation of immunoglobulins and other serum proteins by immunological methods Clin. Chem. Acta 1968 Sep, 22 (1):43.
4. Milford Ward A. (1990). Protein reference Unit Handbook of clinical Immunochemistry. 5th edition, Publ.PRU Publications, Sheffield, UK.
5. Alfonso E. (1964) Quantitative immunoelectrophoresis of serum proteins Clin. Chem. Acta 10, 114-122.

Legend of the symbols used on the labels:/

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AMS S.p.A. – Sede legale e stabilimento