# ENZYMATIC HbA1c CALIBRATOR

Calibrator for Glycated Hemoglobin HbA1c assay

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IVD

## INTENDED USE

This product is intended for the calibration of Enzymatic HbA1c assay.

## COMPOSITION

**2 Levels (lyophilized) x 0.5 ml**

Stabilized hemolysates prepared from human whole blood.

## Reference values

Calibrator values are assigned with IFCC standardized method and NGSP referenced materials. Conversion between DCM (Designated Comparison Method) are provided with Master Equation. Values are reported in the attached Value Sheet, lot dependant.

## AUXILIARY EQUIPMENT

- Enzymatic HbA1c REF GD5418 00
- 0.5 ml pipette
- Deionised water

## PREPARATION

Reconstitute each calibrator using 0.5 ml of deionised water. To ensure complete reconstitution, equilibrate vials at room temperature for 30 minutes before first use.

## STORAGE AND STABILITY

Store at 2-8 °C. The calibrators are stable until expiration date if stored at the indicated conditions, protected from light and heat.

## After reconstitution:

The reconstituted calibrators are stable 15 days at 2-8 °C. It is suggested to lyse all the material by following procedure described in Enzymatic HbA1c (REF GD5418 xx) and store at -20 °C the hemolysate. Avoid freeze-thawing cycle.

## PROCEDURE

The lyophilized HbA1c calibrators should be treated in the same manner as specimens, hemolysate procedure included. Follow the instructions for use accompanying the kits and the applications for the instruments.

## Test limitation

Improper pipetting, inadequate mixing and erroneous instruments calibration may cause inaccurate results.

## PRECAUTIONS IN USE

All reagents contain inactive components such as preservatives (Sodium azide or others), surfactants etc. The total concentration of these components is lower than the limits reported by 67/548/ECC and 88/379/EEC directives about classification, packaging and labelling of dangerous substances. However, the reagents should be handled with caution, avoiding swallowing and contact with skin, eyes and mucous membranes.

The use of the laboratory reagents according to good laboratory practice is recommended.

The reagents from human donors have given negative results to anti-HIV 1/2, HBsAg and anti-HCV. It is recommended to handle with caution.

## Waste Management

Please refer to local legal requirements.

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