GLICOTEST HbA1c Calibrator/Calibratore

REF

GD5426 00

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IVD

INTENDED USE

This product is for the purpose of calibrating results in the quantitative determination of human hemoglobin A1c (HbA1c) in blood by automated immunoassay.

COMPOSITION

4 Levels (lyophilized) x 0.5 ml

Stabilized hemolysates prepared from human erythrocytes.

Reference values

The % HbA1c value of the calibrators were obtained by assaying representative samples of the entire lot against materials referenced to NGSP values using the GLICOTEST HbA1c REF GD5424 00 kit. The values are stated on each calibrator vial label.

AUXILIARY EQUIPMENT

- Glicotest HbA1c REF GD5424 00
- 0.5 ml pipette
- Deionized water

PREPARATION

Reconstitute each calibrator using 0.5 ml of deionized water. Gently mix for 10 minutes, or until all material has dissolved. Subdivide the reconstituted calibrator into aliquots of 0.1 ml and seal tightly.

STORAGE AND STABILITY

Store at 2- 8°C. The calibrators are stable until expiration date if stored at the indicated conditions, protected from light and heat.

After reconstitution:

The reconstituted calibrators are stable 30 days at 2-8 °C and 3 months at -20 °C. Avoid freeze-thawing cycle.

PROCEDURE

The lyophilized HbA1c calibrators should be treated in the same manner as specimens, hemolysate procedure included.

Follow the instructions for use accompanying the kits and the applications for the instruments.

Test limitation

Improper pipetting, inadequate mixing and erroneous instruments calibration may cause inaccurate results.

PRECAUTIONS IN USE

All reagents contain inactive components such as preservatives (Sodium azide or others), surfactants etc. The total concentration of these components is lower than the limits reported by 67/548/ECC and 88/379/EEC directives about classification, packaging and labelling of dangerous substances. However, the reagents should be handled with caution, avoiding swallowing and contact with skin, eyes and mucous membranes.

The use of the laboratory reagents according to good laboratory practice is recommended.

The reagents from human donors have given negative results to anti-HIV 1/2, HBsAg and anti-HCV. It is recommended to handle with caution.

Waste Management

Please refer to local legal requirements.

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