

Glycated Haemoglobin A_{1c}

Chromatographic-spectrophotometric
Ion exchange method

Code HB031 20 tests

Stored at 2-25°C.



Clinical Significance

Diabetes Mellitus is a chronic disease characterized by a hyperglycemia. The consequences are metabolism disorders of carbohydrates, lipids and proteins. The risk of complications associated with diabetes, including nephropathy, retinopathy and cardiovascular diseases, increases in patients with poor metabolic control. In the diabetic patients, where blood glucose levels are elevated, HbA_{1c} is formed as a consequence of the non-enzymatic glycation of the N-terminus of the β -chain of haemoglobin molecule. The level of HbA_{1c} is proportional to the level of glucose in the blood and has been widely accepted as an indicator of the mean daily blood glucose concentration over the preceding 6-8 weeks. It is therefore, a long-term indicator of diabetic control, whereas, the measurement of blood glucose is only a short-term indicator.

Principle

After preparing the hemolysate, where the labile fraction is eliminated, hemoglobins are retained by a cationic exchange resin. Hemoglobin A_{1c} is specifically eluted after washing away the HbA_{1a+b} fraction, and is quantified by direct photometric reading at 415 nm.

Reagents

Reagent 1 (1x6 ml)	Potassium biphthalate Detergent, pH 5.0	50 mmol/l
Reagent 2 (1x 47 ml)	Phosphate buffer, pH 6.5 Sodium azide	48 mmol/l 0.95 g/l
Reagent 3 (2x180 ml)	Phosphate buffer, pH 6.4 Sodium azide	72 mmol/l 0.95 g/l
Reagent 4 20 pcs	Microcolumns: Contain a pre-weighed amount of resin equilibrated with phosphate buffer.	

Use only columns and reagents 2 or 3 of the same lot number. Because the reagents contain sodium azide, it is advisable not to pipet them by mouth.

For *in vitro* diagnostic use only.

Preparation

Reagents are ready to use.

Bring the column and reagents to room temperature (21-26°C), a few minutes before use.

Column preparation: Always remove the upper cap from the column first and the lower cap second. Then, using the rounded end of a pipette, push the upper disc down to the resin surface taking care not to compress it. Let the column drain completely to waste.

Storage and stability

All the components of the kit are stable at 2-25°C up to the date of expiration as specified.

Additional equipment

- Spectrophotometer or colorimeter measuring at 415 nm
- Matched cuvettes 1,0 cm light path
- Test tubes (16x160 mm)
- General laboratory equipment

Samples

Whole blood collected in heparin or EDTA. Stable for at least 7 days at 2-8°C

Procedure

1. Wavelength 415 nm (405-425); Temperature 21-26°C; Cuvette 1 cm light path.
2. Adjust the instrument to zero with distilled water.
3. **Hemolysate preparation and labile fraction elimination**

Pipette into a test tube:

Blood	50 μ l
Reagent 1	200 μ l

Shake thoroughly and let it stand at room temperature for 10-15 minutes. The hemolysate will be used in steps 5 and 6.

4. Prepare a column: see preparation
5. **Separation and reading of HbA_{1c} fraction:**
Carefully pipette on the upper filter:

Hemolysate	50 μ l	Let the column drain to waste
In order to drain any sample residue left above the upper disc, pipette:		
Reagent 2	200 μ l	Let the column drain to waste
Pipette:		
Reagent 2	2.0 ml	Let the column drain to waste
Place the column over a test tube (16 x 100 mm) and add:		
Reagent 3	4.0 ml	Collect the eluate (HbA _{1c} fraction)
Shake thoroughly and read the absorbance (Abs) of the HbA _{1c} fraction at 415 nm against distilled water (Abs _{HbA1c})		

6. **Reading of Hb_{TOTAL}**

Pipette into a test tube (16x160 mm):

Reagent 3	12.0 ml
Hemolysate	50 μ l
Shake thoroughly and read the absorbance (Abs) at 415 nm against distilled water (Abs _{HbTOTAL})	

Calculation

%HbA_{1c}

$$= \text{Abs}_{\text{HbA1c}} / (3 \times \text{Abs}_{\text{HbTOTAL}}) \times 100$$

Reference values

Normal range:	4.2 – 6.2 %
Good controlled diabetic:	6.7 – 9.1 %
Uncontrolled diabetic	> 9.1 %

This range is given for orientation purpose only. Each laboratory should establish its own reference range.

Performance characteristics

Measuring range: from 4.3% (detection limit) to at least 17.0%.

Precision:

	Intra-assay (n=25)		Inter-assay (n=25)	
Mean (mg/dl)	7,2	9,9	7,2	9,9
CV (%)	5,4	6,3	7,3	5,9

Interferences

Erroneous values might be obtained from samples with abnormally elevated quantities of other hemoglobins as a result of either their simultaneous elution with HbA_{1c} (HbF) of differences in their glycation from that of HbA (HbS).

Notes

1. The obtained values are temperature-independent when working in the recommended interval (21 – 26°C). If working temperature is out of range, multiply the obtained value by the corresponding factor showed in the following table:
Working temp. Factor
18°C 1.29
19 °C 1.18
20°C 1.12
27-30°C 0.91
2. The long-term storage of the columns leads to an excessive packing of the resin diminishing the flow rate and lengthening the elution step. To regain the flow efficiency it is advisable - 10 min. before starting the test - to invert the columns to resuspend the contents, place them back to their upright position and let the resin settle for a few minutes
3. Some air bubbles may occasionally appear inside the resin bed. Their presence do not alter the test performance.

Bibliography

1. Bisse E, Abraham EC. J Chromatolog 1985; 344:81-91
2. Tietz NW. Clinical guide to laboratory tests, 3rd ed. Saunders Co., 1999
3. N. Engl. J. Med. 1993; 329: 977-986
4. Young DS. Effects of drugs on clinical laboratory tests, 4th ed. AACCPress, 1995
5. Roberts WL et al. Clin. Chem. 2002; 48: 383-385
3. Goldstain DE et al. Clin Chem 1986; 32: B64-B70.
6. Bry L. et al. Clin Chem 2001; 47:153-163

Langdorp, 12. 2006.