LIPOPROTEINS Controls

Multiparametric controls for lipoproteins assays (human matrix)

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IVD

INTENDED USE

Lipoproteins Controls is a set of multi-analyte controls to be used for monitoring the accuracy and precision of HDL Cholesterol, LDL Cholesterol, Total Cholesterol and Triglycerides assays, with manual and automated procedures.

COMPOSITION

CONTROLS: 3 levels x 3 ml (lyophilized)

Lyophilized human serum base
Preservatives and stabilizers

The concentration of the components are indicated in the attached value sheet and are lot-specific.

PREPARATION

- Open the vial carefully, avoiding any loss of the material and reconstitute with 3 ml of distilled water.
- Replace the rubber stopper, close the vial and leave to stand for 30 minutes before use.
- Ensure that all traces of dry material are dissolved by swirling gently.
- Replace cap immediately after use.

STORAGE AND STABILITY

Controls are stable up to the expiration date stated on the label, if stored unopened at 2-8 °C.

After reconstitution

Once reconstituted, the components of the serum are stable for 7 days at 2-8 °C and 4 weeks at -20 °C, when frozen once. Values may drop up to 10% for LDL Cholesterol when stored for 4 weeks at -20 °C.

PRECAUTIONS IN USE

All reagents contain inactive components such as preservatives (Sodium azide or others), surfactants etc. The total concentration of these components is lower than the limits reported by 67/548/ECC and 88/379/EEC directives about classification, packaging and labelling of dangerous substances. However, the reagents should be handled with caution, avoiding swallowing and contact with skin, eyes and mucous membranes.

The use of the laboratory reagents according to good laboratory practice is recommended.

The reagents from human donors have given negative results to anti-HIV 1/2, HBsAg and anti-HCV. It is recommended to handle with caution.

Waste Management

Please refer to local legal requirements.

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